

Guideline for LC-MS data analysis for facility users

Liquid chromatography–mass spectrometry (**LC-MS**, or alternatively **HPLC-MS**) is an analytical tool that combines the physical separation capabilities of liquid chromatography (or HPLC) with the mass analysis capabilities of mass spectrometry (MS). LC-MS is a powerful technique that has very high sensitivity and selectivity and so is useful in many applications.

Experiment

1. **LC-ESI/MS (LC-MS)2ChScan****: It is a two channel fullscan data acquisition experiment with polarity ES positive and ES negative at cone voltage 30, simultaneously.

Objective: Determination of molecular weights of compounds/analytes separated by liquid chromatography in LC-MS system.

2. **LC-ESI/MS (LC-MS)4ChScan****: It is a four channel fullscan data acquisition experiment with polarity ES positive and ES negative at cone voltage 30 and 60, simultaneously.

Objective: Determination of molecular weights of compounds/analytes separated by liquid chromatography in LC-MS system with In-source fragmentation.

Key Aspects of ESI In-Source Fragmentation (4ChScan)

Structural information: It can provide characteristic fragments for identification, often used when MS/MS is not available. It serves as an alternative to tandem MS (MS/MS) for structural elucidation.

https://www.saiflcknow.org/download_file

SAMPLE COPY OF LC-MS DATA/RESULTS-A

- **SAMPLE COPY OF LC-MS DATA/RESULTS-A**

A sample copy of LC-MS data (2ChScan) is for demo purposes; users can see it. Note: LC-MS data analysis is part of the end-user; it will not be provided by SAIF.

2025-02-15 00:00:00 | [Download](#)

- **SAMPLE COPY OF LC-MS DATA/RESULTS-B**



A sample copy of LC-MS data (4ChScan) is for demo purposes; users can see it. Note: LC-MS data analysis is part of the end-user; it will not be provided by SAIF.

2025-02-15 00:00:00 | [Download](#)

Note:


- The results of the aforementioned LC-MS experiment will be a processed data file in PDF format, and only manual data analysis is possible.
- Data analysis is part of the end-user; it will not be provided by SAIF.

Data interpretation/analysis: LC-ESI (+) MS mode




Data interpretation/analysis: LC-ESI (-) MS mode

Openlynx Report SAIF, CSIR-CDRI, Lucknow
 Sample: 7
 File:19IDKMAUG07
 Description:BETASIL C8 250 X4.6, 5um
 Vial:1:57
 Date:07-Aug-2019
 ID:DKMR60
 Time:22:29:44
 Page 2
 Printed: Thu Aug 08 09:24:04 2019



Openlynx Report SAIF, CSIR-CDRI, Lucknow
 Sample: 7
 File:19IDKMAUG07
 Description:BETASIL C8 250 X4.6, 5um
 Vial:1:57
 Date:07-Aug-2019
 ID:DKMR60
 Time:22:29:44
 Page 32
 Printed: Thu Aug 08 09:24:04 2019



Sophisticated Analytical Instrument Facility & Research. CSIR-Central Drug Research Institute, Sector 10, Jankipuram Extension, Sitapur Road, Lucknow 226031, Uttar Pradesh, INDIA

Electrospray ionization (ESI) source adduct formation occurs when analyte molecules bind with ions (e.g. H^+ , Na^+ , NH_4^+) in the spray plume, creating adduct ions like $[M+H]^+$, $[M+Na]^+$, $[M+NH_4]^+$ or $[M-H]^-$, $[M+Cl]^-$, $[M+CH_3COO]^-$. These adducts are crucial for determination of molecular weights of analytes. User may use <https://tmsdatabase.cdri.res.in/AductionCalculator> for possible adduct formation and make a table for as shown below.

Example-1

Table -1: Chromatographic retention time and detection of compound with adduct formation.

S.No.	RT	M.W.	ESI (+)		ESI (-)		
			$[M+H]^+$	$[2M+H]^+$	$[M-H]^-$	$[2M-H]^-$	$[M+Cl]^-$
1	6.18	492	493	985	-	-	-
2	7.18	624	-	-	623	1247	659

Note-

- Users can anticipate the probable chemical compounds in a sample based on previously known compounds and their molecular weight information in the literature (written papers on similar samples such as plants, microorganisms, and animals).
- The user should do extensive literature search so that the maximum number of compounds can be identified in sample. Without literature support, LC-MS data analysis is not possible.
- After obtaining LC-MS data and molecular weight information, the user need LC-HRMS and LC-HRMS/MS data to determine chemical formula and identify the chemicals detected in MS.
- The user can utilize the ESI/APCI ion calculator to confirm the molecular weight of chemicals observed during the LC-MS run. <https://tmsdatabase.cdri.res.in/AductionCalculator>
- Data analysis is part of the end user and will not be provided by SAIF.

Reference:

- Singh A, Ali M, Mishra DK, Kanojiya S. Comparative Profiling of Withanolides in Tissues of *Withania somnifera* and *Datura metel* Using UHPLC-MS/MS. *Chem Biodivers*. 2025 Dec;22(12):e01319. doi: 10.1002/cbdv.202501319. Epub 2025 Aug 22. PMID: 40845291.
- Singh Y, Nimoriya R, Rawat P, Mishra DK, Kanojiya S. Structural Analysis of Diastereomeric Cardiac Glycosides and Their Genins Using Ultraperformance Liquid Chromatography-Tandem Mass Spectrometry. *J Am Soc Mass Spectrom*. 2021 May 5;32(5):1205-1214. doi: 10.1021/jasms.1c00017. Epub 2021 Apr 5. PMID: 33818079.